



**ABEL TECH**  
SUPPORT YOUR LIFE  
AN A+ COMPANY

Model No.: DSCO-T05

The most compact size camera head design for 720P model

**Handy, Tiny Head & Slim Type USB Dental Camera w/ High Quality & Inexpensive**  
It is also a good tool for Telehealth w/ Intraoral & extra-oral functions. **Competitive Cost Model**

2MP/720P



8 pcs White LEDs



8.5 mm Height    12 mm Width



USB into Android Tablet PC for TeleHealth (TeleMedicine) to talk with Doctor or AI System



Direct to PC



USB to PC for Dental PMS / Image

USB Cable fixed into Camera Body tightly, no more special Video connector; special connector is used to loose and high cost once changed

The capture button working w/ Image Software (Dexis, Carestream, Dental 4 Windows, Patterson Eaglesoft Imaging, VixWin, Schick CDR DICOM, ExaminePro, Praktika, Dentimax, XDR, Tiger View, Sopro, Cliniview, Digora for Windows, Mediasuite, Romexis, Sidexis, Ultimo, Examine Pro via integration software w/ TWAIN Driver

**DEXIS**

Carestream

schick  
by Sirona

PATTERSON  
DENTAL



Iodak | Dental Systems

MiPACS

eagle soft

GENDEX  
Always by your side

PLANMECA

SOPRO  
imaging

TigerView

XDR

DentiMax

mediasuite<sup>sql</sup>  
Leading Digital Imaging Software

CliniView™  
Dental Imaging

DentalEye

APTERYX  
custom software

Umbie DentalCare  
PRACTICE MANAGEMENT MADE EASY

Praktika  
online dental software