

Model No.: DSCO-T41

# 4K VIDEO INTRAORAL CAMERA



◆ Compact Camera Head w/12 LED, 4 white LED, 8 blue LED; under Blue LED mode, it is easier to find Dental Plaque & Calculus, and Decayed teeth



◆ USB Cable fixed into housing, more reliable and stable



◆ Embedded ASCII Code (HID) Control for Buttons & Auto LED OFF/ ON



● USB to PC for Dental PMS / Image

## < Features >

- ◆ Pioneer 4K Video Intraoral Camera w/ USB type C
- ◆ Tiny & Narrow Camera Head w/ 12 Mega Pixel CMOS Sensor
- ◆ Advanced Optical Camera Lens to work with Sony Sensor
- ◆ Following Video Resolution alternative
  - 4K / 3840x2160 pixel with 25 fps
  - 1080P / 1920x1080P pixel with 30 fps
  - 720P / 1280x720P pixel with 60 fps
- ◆ Super Video Image Quality w/ 4K Video Preview
- ◆ Equipped w/ ASCII code (HID) to work w/ special Image Software
- ◆ Designed & Made in Taiwan
- ◆ 110 ° View Angle w/ long focal distance
- ◆ Focus for one tooth, one row teeth, whole mouth, even face
- ◆ Physical buttons: Capture / Zoom IN & OUT / Freeze / LED OFF & ON
- ◆ USB 2.0 w/ long USB cable fixed into camera body