



# G5

INNOVATIVE DESIGN  
WITH EXQUISITE CRAFTSMANSHIP

...

# Main Control Instrument Panel

Fully Displays the Advantages  
of Intelligent Personalized Customization

Brand New ABS Precision Plastic  
Does Not Yellow Over Long-Term Use



## Button Control Panel

Sensitive and durable, passed tens of thousands of button press tests

Icons are intuitive and clear, easy to operate

Clear functional partition, enhances Control Precision

Three sets of seat memory functions, intuitive pressure gauge



## Pure copper phone connector

More durable



## More Details



Stainless steel rotating tray

Easy to take and place, worry-free cleaning

# Secondary Control Panel

Four-hand Efficient Operation  
Provides More Efficient and Convenient  
Functions and Operational Space

- **Button-type Control Panel**

Chair position adjustment button

Mouthwash injection, heating, flushing button

Spitting position, reset, oral light

- **Multi-joint, Flexible Position Adjustment**



# Induction Oral Lamp

Three Major Intelligent Modes  
to Meet Different Treatment Needs

**2** types

Operation Mode  
Sensing/Manual

**3** types

Color Temperature  
White/Yellow/Mixed

**55000**

lux (max)  
Brightness control

**4000**

-6000k  
CT control



# Multifunctional Foot Switch

Noise-free Control

Tested for 300,000 Uses, Durable

## Chair position control

Lift/spit position/reset

## Suction

Mouthwash switch

## Water and air operation

Water spray/air spray

## Debris blow button

Oral light control



**Large pedal**  
Press-type design

**Fixed handle**  
For easy mobility

# Side Cabinet

Can Rotate 60°

Four-hand Operation is More Flexible

- Imported Water and Gas Tubes
- Imported Solenoid Valve
- Independent phone tube
- Water and gas control
- Heater: Warm water supply



# More Details



**Integrated suction filter**  
Effective filtration, easier cleaning




**Ultra-deep rotatable ceramic sputum cup:** can rotate 90°  
easy to disassemble, durable sputum cup filter and cap



**Double bottle water supply system**  
Independent water supply to prevent cross-contamination.

# Ergonomic Design

Human-centered design concept  
Better consideration of patient needs and comfort



Double-joint headrest



Special design backrest flat  
ultra-thin wide wing

# Streamlined exterior design


Pursue ultimate smoothness and simplicity,  
make the diagnosis and treatment environment  
more comfortable and pleasant



# Light luxury cushion

Exquisite and comfortable seat

The first step to a warm and comfortable atmosphere



Breathability ★★★★★


Softness ★★★★★

Durability ★★★★★

Selected microfiber leather with delicate lychee texture



Thickened cushion  
for double comfort



Available in multiple colors